

# **HAEMEDIC<sup>®</sup> Plus**



**Safety lancets for capillary blood sampling**

**HaeMedic**

# HAEMOLANCE<sup>®</sup> Plus

## Safety Lancets

### Features and benefits

- Needle retracts automatically and locks out for extra safety, protecting patient and health professionals from needlestick injuries and prevent cross-infection.
- Ergonomic and easy to handle.
- HAEMOLANCE Plus saves time. It is quick and simple to use and requires little instruction.
- HAEMOLANCE Plus can not be reloaded or reused.
- Simple and safe for both patient and user!

### How to use HAEMOLANCE Plus



1. To remove the protective cap turn and twist in one movement (a full turn)



2. Place lancet against side of fingertip and press button firmly.

	Name	Needle Size	Puncture Depth	Product No.	Pcs/Box	Transport and ordering pack.
	Paediatric	Blade	1,2 mm	53164	150 pcs	24 boxes
	Max Flow	Blade	1,6 mm	53169	150 pcs	24 boxes
	High Flow	18~G	1,8 mm	51361 51362	150 pcs 50 pcs	24 boxes 54 boxes
	Normal Flow	21~G	1,8 mm	53152 53151	150 pcs 50 pcs	24 boxes 54 boxes
	Low Flow	25~G	1,4 mm	53156 53155	150 pcs 50 pcs	24 boxes 54 boxes
	Micro Flow	28G	1,6 mm	53159 53158	150 pcs 50 pcs	24 boxes 54 boxes

#### HaeMedic

Företagaregatan 18,  
266 32 Munka Ljungby  
Sweden  
Phone: +46 431-43 20 30  
Fax: +46 431-43 00 30

#### HaeMedic

P. O. Box 806  
Wayzata, MN 55391-0806  
USA  
Phone: 763 404 9651  
Fax: 763 404 5997

E-mail: [info@haemedic.se](mailto:info@haemedic.se) [www.haemedic.com](http://www.haemedic.com)

Made in Sweden by Arta Plast AB,  
Antennvägen 1A, 135 48 Tyresö, Sweden

CE 0413