

●●●●●●

PROLANCE

Safety lancets for capillary blood sampling



PROLANCE

Safety lancets for capillary blood sampling

- Needle retracts automatically and locks out for extra safety, protecting patient and health professionals from needle stick injuries.
- Six different sizes available allowing health professionals to derive the proper amount of blood for any capillary test.
- Prolance is quick and simple to use, requiring little instruction.
- Prolance cannot be reloaded or reused, protecting both patients and health professionals.
- Double spring action to provide greater patient comfort by firing and retracting the lancet with spring action
- Designed to provide controlled needle depth allowing for greater patient comfort.
- Sterilized by e-beam, assuring sterility for the patient.
- Ergonomic and easy to handle with side ribs allowing for safe handling.







How to use PROLANCE



1. Holding the Prolance by the edges remove protective cap by turning 360 degrees or more.



2. Place the tip of the lancet firmly but gently against the finger approximately 2 mm from nail. Press and release button in one quick movement. Dispose of properly after use.

	Name	Needle Size	Puncture Depth	Product No.	Pcs/Box	Transport and ordering pack.
	Prolance Pediatric	Blade	1.2 mm	7609	200 pcs	18 boxes
	Prolance Max Flow	Blade	1.6 mm	7620	200 pcs	18 boxes
	Prolance High Flow	18≈G	1.8 mm	7601	200 pcs	18 boxes
	Prolance Normal Flow	21≈G	1.8 mm	7594	200 pcs	18 boxes
	Prolance Low Flow	25≈G	1.4 mm	7599	200 pcs	18 boxes
	Prolance Micro Flow	28G	1.6 mm	7604	200 pcs	18 boxes