

**SHIPPING CONFIGURATION OF HAEMOLANCE PROFILE® SAFETY LANCETS**

**Shelf boxes:**

PRODUCT	NEEDLE /BLADE	PENETRATION DEPTH	COLOR CODING	CAT No.	SHELF BOX
Haemolance® profile Low Flow	Needle 21G	2.0 mm	Blue	7510	150 pcs
				7509	50 pcs
Haemolance® profile Normal Flow	Needle 21G	2.0 mm	Green	7501	150 pcs
				7500	50 pcs
Hamolance® profile High FFlow	Needle 21G	2.0 mm	Yellow	7516	150 pcs
				7515	50 pcs

**Transport packaging:**

Shelf boxes are packed in:

- Intermediary box: 31.2 x 40.8 x 45.2 cm – 50 pcs/Box
- Intermediary box: 38.5 x 38.5 x 28.0 cm – 150 pcs/Box
- Transport cases

Shelf box:	150 pcs	50 pcs
Transport case:	30 shelf boxes/ 4,500pcs	80 shelf boxes/ 4000 pcs
Transport case weight:	7 kg	7 kg

**Pallet Configuration:**

Pallet	Transport Cases	Shelf boxes
1 pallet 80 x 120 cm	42	1,260 boxes of 150 pcs/ 96,000 pcs
	24	1,920 boxes of 50 pcs/ 189,000 pcs
1 seapallet High cube (USA)	63	1,890 boxes of 150 pcs/ 283,500 pcs
	40	3,200 boxes of 50 pcs/ 160,000 pcs

**Min quantity of single order:**

- Min order: 1 transport case = 24 boxes of 200 pcs or 48 boxes of 100 pcs (4,800 pcs)
- Regular sale: orders for 1 pallet quantity would be appreciated = 720 boxes of 200 pcs or 1,440 boxes of 100 pcs (144,000 pcs)